
Invacare® Robin® EVO

en Ceiling Hoist

User Manual 2

da Løftemotor

Bruksanvisning 28

no Takmontert løfter

Bruksanvisning 54

sv Taklyft

Bruksanvisning 80

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and
saved for future reference.

Yes, you can.®

Contents

1 General	3
1.1 Introduction	3
1.1.1 Symbols in this Document	3
1.2 Service Life	3
1.2.1 Additional information	3
1.3 Limitation of Liability	3
1.4 Warranty Information	3
1.5 Compliance	3
1.5.1 Product-specific standards	4
2 Safety	5
2.1 General safety information	5
2.2 Safety Information on Accessories	6
2.3 Safety information on electromagnetic interference	6
2.4 Labels and symbols on the product	6
2.4.1 Label location	6
2.4.2 Identification label	7
2.4.3 Read the User Manual Label	7
2.4.4 Symbols on the Lifting Module's Display	7
2.4.5 Symbols on the Hand Control	8
3 Product Overview	9
3.1 Intended Use	9
3.2 Main parts of the ceiling hoist	9
3.3 Options	10
4 Setup	11
4.1 Included items	11
4.2 Before Use	11
4.3 Installing the ceiling hoist	11
4.4 Hand control plug	12
5 Usage	13
5.1 General safety information	13
5.2 General operating information	13
5.3 Operating the ceiling hoist	13
5.4 Display and Indicator Light	13
5.5 Indicator light and audio signal	13
5.6 Charging the Battery	14
5.7 Emergency stop	15
5.8 Emergency lowering	15
5.9 Resetting the Emergency Stop	15
5.10 Resetting the System	15
5.11 Mechanical Emergency Lowering	16
6 Patient Transfer	17
6.1 General Safety Information	17
6.2 Attaching the Sling	17
6.3 Lifting Information	18
6.4 Lifting to and from a seated position	18
6.5 Lifting to and from a lying position	19
7 Maintenance	20
7.1 General maintenance information	20
7.2 Daily inspections	20
7.3 Cleaning and Disinfection	20
7.3.1 General Safety Information	20
7.3.2 Cleaning Intervals	20
7.3.3 Cleaning instructions	20
7.3.4 Disinfection Instructions	20
7.4 Service interval	21
7.4.1 LOLER Statement	21
7.5 Battery Maintenance	21
8 After Use	22
8.1 Transport and Storage	22
8.2 Dismounting the ceiling hoist	22
8.3 Disposal	22
8.4 Reconditioning	22
9 Troubleshooting	23
9.1 Identifying and Resolving Faults	23
10 Technical Data	25
10.1 Dimensions and weight	25
10.2 Electrical System	25
10.3 Environmental conditions	25
11 Electromagnetic compatibility (EMC)	26
11.1 General EMC information	26
11.2 Electromagnetic emission	26
11.3 Electromagnetic Immunity	26
11.4 EMC test specifications	27

1 General

1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Only use this product if you have read and understood this manual. Seek additional advice from a healthcare professional who is familiar with your medical condition and clarify any questions regarding the correct use and necessary adjustment with the healthcare professional.

Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

Previous product versions may not be described in this manual's current revision. If you require assistance, please contact Invacare.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the product, for example product safety notices and product recalls, contact your Invacare distributor. See addresses at the end of this document.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

1.1.1 Symbols in this Document

Symbols and signal words are used in this document and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING!

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION!

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

NOTICE!

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations, and information for efficient, trouble-free use.

Other Symbols

(Not applicable for all manuals)

Triman

Indicates recycling and sorting rules (only relevant for France).

UKRP

Indicates if a product is not manufactured in the UK.

1.2 Service Life

The expected service life of this product is eight years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

1.2.1 Additional information

The expected service life is based on an estimated average of 4 lifting cycles per day.

A lifting cycle consists on the raising and lowering of the patient lift straps, for the same distance, in both directions.

1.3 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

1.4 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

1.5 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

We comply with the current environmental legislations WEEE and RoHS.

1.5.1 Product-specific standards

The product has been tested and conforms to ISO 10535 (Hoists for the transfer of persons with disabilities) and all related standards.

For further information about local standards and regulations, contact your local Invacare representative. See addresses at the end of this document.

2 Safety

2.1 General safety information

This section of the manual contains general safety information about your product. For specific safety information, refer to the appropriate section of the manual and procedures within that section.

WARNING!

Risk of injury or damage

- Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, Invacare provider or qualified technician before attempting to use this product.
- Do not make any unauthorized alterations or modifications to the product.

WARNING!

The maximum safe working load must not be exceeded

- Do not exceed the maximum safe working load of this product or used accessories like slings, spreader bars etc. See documentation or the labeling for the stated maximum safe working load.
- The component with the lowest load limit determines the maximum safe working load of the entire system.

WARNING!

Risk of injury or damage

Improper use of this product may cause injury or damage.

- Do not attempt any transfer without approval of the patient's healthcare professional.
- Read the instructions in this user manual and observe trained personnel performing transfer procedures. Then practice transfers under supervision and with a capable person acting as a patient.
- Special care must be taken with people with disabilities who cannot cooperate while being transferred.
- The ceiling hoist is advised only to be used to lift or move a patient.
- The red emergency strap must be adjusted to the caregiver's reach, and must not be removed.
- If the hoist fails, lower patient onto a suitable surface using the red emergency strap before continuing the transfer by a different method.

WARNING!

Risk of injury or damage

Inappropriate handling of cables can cause electrical shock and product failure.

- Do not kink, shear or otherwise damage the cables of the product.
- Make sure that no cables are jammed or damaged, when the product is used.
- Make sure the wiring is correct and connections are proper.
- Do not use unauthorized equipment.

WARNING!

Risk of injury or damage

Excessive moisture will damage the product and may cause electrical shock.

- The patient lift can be used in a bath or shower area but must NOT be used under the shower. The patient must be transferred to a shower chair or use other means for showering.
- If the patient lift is used in a moisture environment, ensure that the patient lift is wiped clean of any moisture after use.
- Do not charge the product in a moisture environment or with wet hands.
- Do not store the product in a damp area or damp condition.
- See 10.3 *Environmental conditions*, page 25.

WARNING!

Risk of injury or damage

Ignition sources can cause burns or fire.

- Patient transfer must be performed with a safety space between the hoist and possible ignition sources (heater, stove, fireplace, etc.)
- The patient and assistants must not smoke during transfer.
- The sling must not be placed over heat sources (heater, stove, fireplace, etc.)

WARNING!

Risk of injury or damage

To avoid injury or damage when operating the product:

- Close supervision is necessary when the product is used near children or pets.
- Do not let children play with the product.

CAUTION!

Risk of injury or damage

The product can get hot when exposed to sunlight or other heat sources.

- Do not expose the product to direct sunlight for prolonged periods.
- Keep the product away from heat sources.

NOTICE!

Accumulation of lint, dust and other dirt can impair the product.

- Keep the product clean.

! NOTICE!
 The electronics may only be serviced by a qualified technician.
 — Do not remove the top cover.

2.2 Safety Information on Accessories

! CAUTION!
Risk of injury
 Non-original or incorrect accessories may affect the function and safety of this product.
 — Due to regional differences, refer to your local Invacare website or catalog for available accessories or contact your Invacare provider.
 — See the manual delivered with the accessory for further information and instructions.
 — Use only original accessories for the product in use. Under certain circumstances, the use slings from other manufacturers is possible. See additional information in this section.

! CAUTION!
Compatibility of slings with attachment system
 Invacare uses a common attachment system based on hooks and loops. Loops on the slings are attached to hooks on the lift. Therefore, suitable slings of other manufacturers can be used on this lift as well.
 — Use only slings with loop attachments that are suitable for hooks as attachment points.
 — Do not use slings designed for "Keyhole or Clip Attachment Systems" or "Tilting Frame Systems".

To select the appropriate sling, a risk assessment must be performed by a healthcare professional. The risk assessment must consider:

- The patient's weight, size, physical ability and medical condition.
- The type of transfer and the environment.
- The compatibility to the other lifting equipment used.


2.3 Safety information on electromagnetic interference

! WARNING!
Risk of malfunction due to electromagnetic interference
 Electromagnetic interference between this product and other electrical equipment can occur and disturb the electrical adjustment functions of this product. To prevent, reduce or eliminate such electromagnetic interference:
 — Only use original cables, accessories and spare parts, to not increase electromagnetic emission or reduce electromagnetic immunity of this product.
 — Do not use portable radio frequency (RF) communications equipment closer than 30 cm to any part of this product (including cables).
 — Do not use this product near active high-frequency surgical equipment and the RF shielded room of a system for magnetic resonance imaging, where the intensity of electromagnetic disturbances is high.
 — If disturbances occur, increase the distance between this product and the other equipment or switch it off.
 — Refer to the detailed information and follow the guidance in chapter 11 Electromagnetic compatibility (EMC), page 27.

! WARNING!
Risk of malfunction
 Electromagnetic interference may cause improper operation.
 — Do not use this product adjacent to or stacked with other electrical equipment. If such use is necessary this product and the other equipment must be closely observed to verify that they are operating normally.


2.4 Labels and symbols on the product

2.4.1 Label location


(A)	SWL — Maximum Safe Working Load
(B)	Model name

Ⓒ	Identification Label
Ⓓ	Read the User Manual
Ⓔ	HMI Label (only applicable for Denmark)

2.4.2 Identification label

The identification label contains the main product information, including technical data.

Symbols	
	Medical device
	European Conformity
	UK Conformity Assessed
	Manufacturer
	Date of manufacture
	Max. safe working load
	Reference number
	Serial number
	Unique Device Identification
	Class II equipment
	Type BF applied part
	WEEE conform

Abbreviations for technical data:

- I_{in} = Incoming Current
- U_{in} = Incoming Voltage
- Int. = Intermittence
- DC = Direct Current
- AC = Alternating Current
- Max = maximum
- min = minute

For more information about technical data, see *10 Technical Data*, page 25.

2.4.3 Read the User Manual Label

Read the User Manual before using this product and follow all instructions for safety and use.


The background colour of this symbol is blue on product labels.

2.4.4 Symbols on the Lifting Module's Display

	Overload
	Strap loose or jammed inside the lifting module
	Strap length reference lost
	Hand control unplugged or malfunction
	Failure of the internal weight sensor
	Failure of the internal motor sensor
	Failure of the internal strap length sensor
	Over-current Motor Stall
	Lowering the lifting straps
	Raising the lifting straps
	Strap reset process is successful

See *9.1 Identifying and Resolving Faults*, page 23.

2.4.5 Symbols on the Hand Control

		<p>Raise the lifting straps, see 5.3 <i>Operating the ceiling hoist, page 13</i></p>
		<p>Lower the lifting straps, see 5.3 <i>Operating the ceiling hoist, page 13.</i></p>

3 Product Overview

3.1 Intended Use

The ceiling hoist is a battery-powered transfer device and is intended to transfer and position an individual from one resting surface to another.

For example:

- Between a bed and a shower trolley
- From a wheelchair to a bed and vice versa

The maximum safe working load is stated in *10 Technical Data, page 25*.

The ceiling hoist is designed to be used indoors in hospitals, nursing facilities and domestic areas in combination with compatible ceiling tracks or gantry.

The ceiling hoist can be a solution in places with limited floor space.

A healthcare professional or private person who has received proper training is the intended operator of this product.

Intended users

The intended user is a completely or partially immobile person.

Indications


The ceiling hoist is indicated for persons who are completely or partially immobile, and unable to be transferred manually.

All position changes are possible without assistance of the patient.

There are no known contra-indications for this product.

If continuously exposed to high humidity, high condensation and corrosives (e.g. chlorine or ammonia gases) the product may be impaired and a reduction in lifetime must be expected.

3.2 Main parts of the ceiling hoist

Ⓐ	Suspension pin
Ⓑ	Display
Ⓒ	Mechanical emergency lowering plugs (Two — one on each side)
Ⓓ	Red emergency strap
Ⓔ	Indicator light
Ⓕ	Lifting straps
Ⓖ	Strap hooks
Ⓗ	Hand control
Ⓘ	Charging station
Ⓢ	Power supply adapter

3.3 Options


Due to regional differences, refer to your local Invacare website or catalog for details about available options or contact your Invacare provider.

- **Two** and **Four** point spreader bar
- Winnicare **Two** point Spreader Bar with an embedded Medical Grade Scale
- Extension kit (extension straps, extension cable for the hand control and extension strap for the emergency strap)
- Sling models for a “Loop and Coat Hanger” Spreader Bar System:
 - Full body support slings – with / without head support
 - Slings for dress / toileting – with / without head support
 - Amputee Slings
 - Slings for walking training

4 Setup

4.1 Included items

(A)	Ceiling hoist (including two lifting straps with hooks and a red emergency strap)
(B)	Hand control
(C)	Charging station
(D)	User manual
(E)	Power supply adapter

4.2 Before Use

Before first use, all features of the hoist system must be inspected by specialized service technicians.

Perform a visual inspection of the ceiling hoist. If the packaging is damaged upon arrival, the ceiling hoist parts must be thoroughly examined for visible damage or defects. If there is any suspicion of damage, the ceiling hoist must not be used before specialized service technicians have approved it.

The emergency stop is activated during shipment. Reset the emergency stop before first use. Refer to *5.9 Resetting the Emergency Stop, page 15*.

It is possible that the hoist's battery is discharged on account of self-discharge. Charge the hoist's battery for at least 10 hours before first use. Refer to *5.6 Charging the Battery, page 14*.

4.3 Installing the ceiling hoist

The ceiling hoist must be mounted to a rail system to perform its intended use. The rail system must be installed and approved according to ISO 10535 by a qualified technician.

This ceiling hoist requires a trolley, holding it inside the rail system. The trolley must be pre-installed by a qualified technician.

Refer to EC-Track installation manual.

WARNING!


Risk of injury or damage

Non-original or incorrect rail system may affect the function and safety of this product.

- If a rail system from another manufacturer is intended to be used with the Invacare ceiling hoist, a risk assessment for the combination/compatibility of both products shall be performed by the responsible person/organization for the combination.

To mount the ceiling hoist to the trolley, do the following:


1. Locate the hoist trolley (B) at the $\varnothing 32$ mm round recess in the rail (A).

2. Lift the hoist and push the suspension pin through the trolley opening.
3. Turn the hoist by 90 degrees to fix it in the trolley.


Keep pushing the device to avoid the free rotation of the pin when fixing in the trolley.

4. Adjust the length of the emergency stop/lowering cord to within the user's reach.


4.4 Hand control plug

Inserting the hand control plug

1. Make sure the hand control plug **A** is correctly aligned with the slot **B**. Insert the push pull self-locking hand control plug into the slot **B**.

Removing the hand control plug

1. Pull the plug **A** out of its slot.

5 Usage

5.1 General safety information

WARNING!

Risk of injury or damage

Before using the lift with a patient, refer to the following safety information and instructions:

- 2 Safety, page 5
- 6 Patient Transfer, page 17.

5.2 General operating information

CAUTION!

If the lifting straps are guided incorrectly in the hoist (e.g. the straps become twisted), the ceiling hoist will automatically turn off.

- Keep the straps straight to ensure they feed into the hoist correctly. The straps are held taut by the weight of the hooks.
- Keep the straps freely suspended when the unloaded hoist is lifting/lowering.

If a user unintentionally traps the hand in the lifting strap during lifting, the hoist will automatically turn off. The hoist stops immediately when a hand gets in contact with the strap inlet.

- Keep your hands away from the strap inlet during lifting.


CAUTION!

The hand control is protected against splashed water (IPX4). Extended contact with water will damage the hand control.

- Do not immerse the hand control in liquids, even for a short time.
- Avoid any contact with water or other liquids.

5.3 Operating the ceiling hoist

The ceiling hoist is in operation mode as soon as a button on the hand control is pressed.

Before use, detach the hand control from its storage location — the charging station —, and put it back when not in use.

Raising and Lowering the Lifting Straps

1. Press and hold button **A** to raise the lifting straps.
2. Press and hold button **B** to lower the lifting straps.

Release the button to stop raising or lowering.

5.4 Traverse along the track

1. Use the two lifting straps to pull the ceiling hoist along the track.


DO NOT use the red emergency strap or the hand control to traverse the ceiling hoist.

The hand control is magnetic and can be placed on the hooks during a transfer.

5.4 Display and Indicator Light

The lifting module features a display **A** and an indicator light **B** to provide information during normal operation, and if an error is detected or corrected.

Normal operation:

- Display — Switches on when: pushing up/down buttons, or performing a system reset. Automatically switches off again.
- Indicator light — permanent/flashing blue light.

The indicator light automatically switches on when a button on the hand control is pressed and switches off when the system switches into standby mode.

Error mode:

- Display — Switches on and shows a symbol for the detected error.
- Indicator light — flashes red.

Error corrected:

- Display — Error symbol disappears.
- Indicator light — blinks green once.


For more information on error symbols, see 2.4.4 Symbols on the Lifting Module's Display, page 7.

5.5 Indicator light and audio signal




Hoist

Status	Indicator light	Audio signal
Off / Stand by	Off	No
Hoist in use (Up/Down)	Flashing Blue	No
Low battery	Flashing Red	Yes


The indicator light is on the underside of the hoist.

Hand Control Battery Level

Charge Level	Indicator light	Functions
50% — 100%	All Three Leds green 	Fully Functional
25% — 50%	First Two Leds green 	Fully Functional
15% — 25%	First Led orange 	Fully Functional
< 15%	First Led blinks orange 	Limited to Strap downward motion
< 10% (Error — battery depleted)	n/a	No functions available

Charge Level	Indicator light	Functions
Charging	0 to 25% — First Led blinks orange 25 to 50% — First Two Leds blink green 50 to 100% — All Three Leds blink green Charging topped — All Three Leds green	No functions available

5.6 Charging the Battery

If the charge status becomes low, an audible indicator sounds and the indicator light on the bottom of the ceiling hoist flashes red (see also 5.5 Indicator light and audio signal, page 13

When this happens, there is usually sufficient power to perform at least one transfer.

 It is not possible to use the hoist while the battery is charging.

1. Finish the current lift in progress.
2. Move the ceiling hoist to the charging station.
3. Assure that the hand control and charging station are clean and dry before charging starts.
4. Place the hand control into the charging station and make sure that the power supply adapter is plugged in. A blinking orange led on the hand control indicates it is charging.

 During the charging process, the LEDs of the hand control are permanently on and blinking.

5. When fully charged, the indicator lights on the hand control turn green (static light). As long as the ceiling hoist is not in use, leave the hand control in the charging station.

 The charging LEDs do not go into sleep mode until the charge has reached full.

6. To use the ceiling hoist, remove the hand control from the charger.


 Keeping the batteries charged guarantees hoist functionality and maintains the battery to ensure a long lifetime.

If the ceiling hoist will not be used for more than 1 week, Invacare recommends disconnecting the charger and pulling the emergency stop.

The charger is connected to a power supply adapter which, in turn, connects to the mains with a mains plug. Make sure that the mains plug is accessible and can be unplugged if required.

5.7 Emergency stop

The emergency stop function is used to stop the patient from lowering or lifting if the ceiling hoist does not stop or react to the hand control.


1. Pull and release the red emergency strap (A) underside of the hoist (B) to activate the emergency stop.

The ceiling hoist stops any lifting or lowering movement immediately.

NOTICE!
When the red emergency strap has been pulled, the hoist functions are deactivated. All visual indications are turned off.
— To activate the hoist functions again, reset the emergency stop. Refer to 5.9 *Resetting the Emergency Stop*, page 15.

5.8 Emergency lowering

The emergency lowering function is used to lower the patient only if the hoist has stopped working.

1. Pull and hold the red emergency strap (A) on underside of the hoist (B) to lower the patient.

 There is an intentional delay of 3 seconds before lowering is initiated.


Keep on pulling, until the patient is lowered to a safe position.

2. Release the red emergency strap to stop lowering the patient.

NOTICE!
When the red emergency strap has been pulled, the hoist functions are deactivated and the blue indicator light turns off.
— To activate the hoist functions again, reset the emergency stop. See 5.9 *Resetting the Emergency Stop*, page 15.

5.9 Resetting the Emergency Stop

WARNING!
— If the emergency stop or emergency lowering must be used, there is a problem with the system. Contact your dealer or Invacare representative for service before using the system again.


When the red emergency strap (B) has been pulled, the reset button (A) pops out so that a red ring is visible. The hoist will not function until the emergency stop or lowering feature is reset.

1. Press the reset button (A) to reset the emergency stop or lowering feature.
2. Perform a System Reset if required, see 5.10 *Resetting the System*, page 15.
 Reset procedure is ONLY needed when the Emergency Lowering is used.
3. Proceed with normal usage.

5.10 Resetting the System

CAUTION!
Risk of Damage to Property
— No weight must be on the device when performing the strap reset.

1. Press and hold down, simultaneously, buttons A and B for three seconds to initiate the system reset and hold the buttons pressed until the following sequence is complete:

a. After three seconds the lifting straps will start to move all the way out (to their maximum length) and then recoil back inside the lifting module.

LCD Display	Indicator light
	Flashing in light blue

b. The lifting straps will stop when they have reached their minimum length.

LCD Display	Indicator light
	Stops flashing in light blue  Indicator light flashes green once when reset is successful. If the procedure is not completed correctly, repeat it.

5.11 Mechanical Emergency Lowering

The Mechanical Emergency Lowering is to be used **EXCLUSIVELY** when one or both motors are not operational.


CAUTION!

A **Mechanical Emergency Lowering** carried out improperly can damage the device.
— Follow the below described procedure carefully.

Before you initiate a **Mechanical Emergency Lowering** try the **Electrical Emergency Lowering**, then do a quick battery charge and retry. Only after this troubleshooting procedures you should initiate a **Mechanical Emergency Lowering**.

Ideally, and if possible, the Mechanical emergency lowering should be done over a bed. If a height-adjustable bed is being used, it can be raised to support the patient during the mechanical emergency lowering procedure.

1. Remove both hole caps A.
2. Insert a 6 mm Allen keys B into the gear motor holes.
3. Rotate both Allen keys B clockwise — simultaneously — until the patient is at a safe height.

WARNING!

If one of the Allen keys on one side is rotated too much, patient can fall.
— If you are unable to rotate both Allen keys simultaneously, alternate rotation on each side.
— Do a maximum of **TEN** full turns on each side.
— Ensure to monitor patient — at all times — during the procedure.

6 Patient Transfer

6.1 General Safety Information

WARNING!

Risk of Injury or Damage

Improper use of this product may cause injury or damage.

- Before transferring to a stationary object (wheelchair, bed etc.) check that its weight capacity can withstand the patient's weight.
- If applicable, the brakes of the stationary object (wheelchair, bed etc.) must be applied before lowering the patient onto or lifting the patient off.
- Always plan the transfer and avoid leaving the patient in the sling unattended.
- Do not lift the patient higher than necessary to perform the transfer.
- Always make sure that the patient will not collide with any stationary objects (wall, furniture etc.) during the transfer.

WARNING!

Risk of Injury or Damage

Objects attached or hung on the device can fall onto the patient and assistant or cause damage due to overload.

- Do not attach or hang any objects on the device.

WARNING!

Risk of Injury or Damage

Damage to parts of the lift caused by impact with the floor, walls or other stationary objects may cause damage to the product and lead to injury.

- Do not allow parts of the lift to impact the floor, walls or other stationary objects.
- Always store the hand control properly when not in use.

WARNING!

Risk of injury

There is a risk of injury to the patient and assistant caused by dangling hooks or the hand control.

- Always be aware of the position of the hooks and hand control during lifting procedures.

WARNING!

Risk of Entrapment or Strangulation

Items in the patient's surroundings can cause entrapment strangulation during lifting. To avoid entrapment or strangulation:

- Before lifting, check that the patient is completely free of his/her surroundings.
- Do not allow the hand control cord or the lifting straps to become entangled around the patient and assistant.

WARNING!

Risk of Entrapment

There is a risk of entrapment between the hooks and the sling.

- Use caution when lifting.
- Never put hands or fingers on or near the hooks when lifting.
- Ensure the patient's hands and fingers are away from the hooks before lifting.

NOTICE!

All transfer procedures described in the following can be performed by one (1) assistant. However, Invacare recommends that the procedures are performed by two (2) assistants whenever possible.

6.2 Attaching the Sling

WARNING!

Risk of Injury

Usage of wrong or damaged slings can cause the patient to fall or cause injury to assistants.

- Use an Invacare approved sling that is recommended by the individual's doctor, nurse or medical assistant for the comfort and safety of the individual being lifted.
- Invacare slings and patient lift accessories are specifically designed to be used in conjunction with Invacare patient lifts.
- After each laundering (in accordance with instructions on the sling), inspect sling(s) for wear, tear, and loose stitching.
- Bleached, torn, cut, frayed, or damaged slings are unsafe and could result in injury. Discard immediately.
- DO NOT alter slings.

WARNING!

Risk of Injury

Improperly attached or adjusted slings can cause the patient to fall or cause injury to assistants.


- Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode).
- DO NOT use any kind of plastic back incontinence pad or seating cushion between patient and sling material that may cause the patient to slide out of the sling during transfer.
- Position the patient in the sling as directed by the instructions provided with the sling.
- Adjustments for patient safety and comfort should be made before moving the patient.

The information in this section is only a general guideline. Refer to the sling user manual for more details about the slings.

A sling with four to six straps designed for mounting on hooks should be used when using the ceiling hoist.


The straps of the sling may be equipped with colour coded loops providing different lengths to place the patient into different positions. Shorter straps at the shoulders will produce a more vertical lift, which will assist with positioning into a chair or wheelchair. By lengthening the straps at the shoulders, a more reclined position can be achieved which is more suitable for a chair to bed transfer. Always match the loop colours of the corresponding straps on each side of the sling for an even lift of the patient.

The sling straps are to be hooked in pairs (one from the back and one from the leg) on their own suspension hook.

1. Place the desired loop of strap on the hook.
2. Pull down the strap until the loop is fully seated at the bottom of hook.
3. Repeat the steps for each of the remaining straps of the sling.

 When attaching the sling to the hooks, and the straps are twisted, the straps will rotate and stretch out, thus ensuring trouble-free operation.

6.3 Lifting Information


WARNING!

Only personnel who have received training regarding the use of lifting equipment and fitting of slings should use the hoist.

- Plan the move. Avoid leaving the client in the sling unattended.
- The ceiling hoist lifts quickly. Before lifting, check that the client is completely free of his/her surroundings.
- The patient's body parts must not be in danger of becoming trapped.
- Check that the red emergency strap, hand control and hand control cable are free of the lifting straps, client and other objects before the ceiling hoist is activated and moved up or down.
- When the ceiling hoist is used correctly, the client should only be lifted to the extent that she/he is clear of the surface and should be moved at this height.


6.4 Lifting to and from a seated position

Follow this procedure when lifting a patient from a seated position (wheelchair, commode, etc.).

1. Position the patient in the sling. Refer to the sling user manual.
2. Traverse the ceiling hoist along the track towards the patient to be lifted.
3. Position the suspension hooks at the same height as the patient's chest and no closer than mid-thigh.
4. Position the suspension hooks so they are parallel to the patient's shoulders.
5. Attach the sling to the suspension hooks. Refer to 6.2 *Attaching the Sling*, page 17.

6.5 Lifting to and from a lying position

Follow this procedure when lifting a patient from a lying position.

1. Apply the sling on the patient. Refer to the sling user manual.
2. Traverse the ceiling hoist along the track towards the patient to be lifted.
3. Position the suspension hooks above the patient to be lifted.
4. Align the suspension hooks with the patient's shoulders.
5. Attach the sling to the suspension hooks. Refer to 6.2 *Attaching the Sling*, page 17.

7 Maintenance

7.1 General maintenance information

CAUTION!

Risk of injury or damage

- Do not perform any maintenance or service procedures while the product is in use.

Follow the maintenance procedures described in this manual to keep your product in continuous service.

Further maintenance and inspection procedures, to be performed by a qualified technician, are described in the service manual of this product. Service manuals can be obtained from Invacare.

7.2 Daily inspections

WARNING!

Risk of injury or damage

Damaged or worn parts can affect the safety of the lift.

- The patient lift should be checked each time it is used.
- Do not use the lift if damage is found or you question the safety of any part of the lift. Contact your Invacare provider immediately and ensure the lift is not used until repairs are made.

Daily Inspection Checklist

- Visually inspect slings, spreader bar (optional), lifting straps and strap hooks. Check all parts for damage or wear.
- Check the emergency stop and emergency lowering function. Check all parts for external damage or wear.
- Check that all hardware and attachment points for damage or wear. Check all parts for external damage or wear.
- Verify that the hand control is functional.
- Check the battery charger.
- Charge the battery every day the lift is used.

7.3 Cleaning and Disinfection

7.3.1 General Safety Information

CAUTION!

Risk of Contamination

- Take precautions for yourself and use appropriate protective equipment.

CAUTION!

Risk of Electric Shock and Product Damage

- Switch off the device and disconnect from mains, if applicable.
- When cleaning electronic components consider their protection class regarding water ingress.
- Make sure that no water splashes to the plug or the wall outlet.
- Do not touch the power socket with wet hands.

NOTICE!

Wrong fluids or methods can harm or damage the product.

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
- Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
- Always make sure that the product is completely dried before taking into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

7.3.2 Cleaning Intervals

NOTICE!

Regular cleaning and disinfection enhance smooth operation, increases the service life and prevents contamination.

Clean and disinfect the product:

- regularly while in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

7.3.3 Cleaning instructions

NOTICE!

- The product does not tolerate cleaning in automatic washing systems, with high-pressure cleaning equipment or steam.

Cleaning the lift

Method: Wipe off with a damp cloth or soft brush.

Solvent/chemicals: Ordinary household cleaning agent and water.

Drying: Wipe dry with a soft cloth.

Cleaning the sling

Refer to the washing instructions on the sling and to the sling manual for cleaning details.

7.3.4 Disinfection Instructions

In Domestic Care

- Disinfectant: We recommend using an alcohol-based surface disinfectant (with 70-90% alcohol).

Read the instructions on the label of your disinfectant. It provides information on the activity spectrum (bacteria, fungi and/or viruses), material compatibility and the correct exposure time.

1. Ensure the surfaces are cleaned before disinfection.
2. Moisten a soft cloth and wipe-disinfect all accessible surfaces and keep them moistened for the exposure time stated on the label of the disinfectant.

3. Allow the product to air-dry.

In Institutional Care

Follow your in-house disinfection procedures and only use the disinfectants and methods specified therein.

7.4 Service interval

WARNING!

Risk of injury or damage

Service must be performed only by a qualified technician.

- Contact your Invacare provider for service.

Service must be performed at least every 12 months unless otherwise stated in local requirements.

The frequency of inspection must be increased, if the product is continuously exposed to high humidity, high condensation and corrosives (e.g. chlorine or ammonia gases), to prevent product impairment.

7.4.1 LOLER Statement

The UK Health and Safety Executive's Lifting Operations and Lifting Equipment Regulations 1998, require any equipment that is used in the workplace to lift a load be subject to safety inspection on a six monthly basis. Please refer to the HSE web site for guidance www.hse.gov.uk.

The person responsible for the equipment must ensure adherence to LOLER regulations.

7.5 Battery Maintenance

WARNING!

Risk of Injury or Damage

Improper handling of lithium batteries may cause injury or damage.

- Do not open the battery housing as damaging the cell or circuitry may develop excessive heat.
- Lithium batteries that are defective, have been damaged or might produce excessive heat or fire, are not allowed for transportation.
- In case the battery turns hot, disconnect it and avoid contact.
- Do not crush or puncture the batteries as it may ignite.

Partial discharge

Lithium-ion batteries perform better when you practice shallow discharges (between 20% to 80%) instead of full discharges. Frequent full discharges can stress the battery and lead to capacity loss over time.

Storage Recommendations

If you plan to store your product for an extended period, charge the battery at least 50% before storage and pull the emergency stop button. Store the product in a cool, dry place within the specified temperature range, and check the battery's charge level periodically during storage. For prolonged storage, ensure the battery retains a charge between 20% to 80% to prevent degradation.

8 After Use

8.1 Transport and Storage


Invacare recommends that the ceiling hoist is always transported and stored in the original packaging. Place the ceiling hoist on a soft surface, such as a cloth or foam rubber mat.

For long term storage, the emergency stop must be activated. This will reduce the discharge of the battery. Refer to 5.7 *Emergency stop, page 15*.

For transport and storage conditions of the ceiling hoist, refer to 10.3 *Environmental conditions, page 25*.

8.2 Dismounting the ceiling hoist

The ceiling hoist can be dismounted for transport or storage.

1. Lift the hoist slightly and push the suspension pin out of its holder in the trolley.
2. Turn the hoist by 90 degrees and pull it out of the trolley.

8.3 Disposal

WARNING!

Environmental Hazard

Device contains batteries.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- DO NOT dispose of batteries in normal household waste.
- DO NOT throw batteries into a fire.
- Batteries MUST be taken to a proper disposal site. The return is required by law and free of charge.
- Only dispose of discharged batteries.
- Cover terminals of lithium batteries prior to disposal.
- For information on the battery type see battery label or chapter 10 *Technical Data, page 25*.

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

8.4 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection
- Cleaning and disinfection

For detailed information, see 7 *Maintenance, page 20* and the service manual for this product.

Make sure that the user manual is handed over with the product.

If any damage or malfunction is detected, do not reuse the product.

9 Troubleshooting


9.1 Identifying and Resolving Faults

Malfunction	Possible Cause	Solution
Audible signal	Overload detected and/or low battery level.	Reduce the load and try again or charge the unit.
Symbol on display and indicator light flashing red. 	Overload detected.	Reduce the load and try again.
Symbol on display and indicator light flashing red. 	Strap is loose or jammed inside the lifting module.	Pull the lifting straps down to stretch them and try again.
Symbol on display and indicator light flashing red. 	Strap length reference lost (after using the emergency lowering function).	Perform a system reset. Refer to 5.10 <i>Resetting the System</i> , page 15.
Symbol on display and indicator light flashing red. 	Hand control is not properly connected or is defective.	Check the connection of the hand control. If a proper connection of the hand control does not solve the problem, contact your Invacare provider for service.
Symbol on display and indicator light flashing red. 	Failure of the internal weight sensor.	Do not use the ceiling hoist and contact your Invacare provider for service.
Symbol on display and indicator light flashing red. 	Strap blockage detected or failure of the internal motor sensor.	Remove any possible blockage at the strap inlet (e.g. twisted straps or other objects touching) and perform a system reset. Refer to 5.10 <i>Resetting the System</i> , page 15. If the error persists, do not use the ceiling hoist and contact your Invacare provider for service.
Symbol on display and indicator light flashing red. 	Failure of the internal strap length sensor.	Do not use the ceiling hoist and contact your Invacare provider for service.
Symbol on display and indicator light flashing red. 	Over-current Motor Stall.	Do not use the ceiling hoist and contact your Invacare provider for service.
Symbol on display. 	Low battery.	Charge the battery. Refer to 5.6 <i>Charging the Battery</i> , page 14.
Ceiling hoist does not respond to the hand control's buttons and the display is blank.	Emergency stop activated.	Check that the emergency stop is not activated. Refer to 5.7 <i>Emergency stop</i> , page 15.
	System does not have power or battery is discharged.	Charge the battery. Refer to 5.6 <i>Charging the Battery</i> , page 14.

Malfunction	Possible Cause	Solution
Lifting straps will only move down and not up.	Hand control defective.	Do not use the ceiling hoist and contact your Invacare provider for service.
	Strap length reference was lost.	Perform a system reset. Refer to 5.10 <i>Resetting the System, page 15.</i>
Lifting straps will only move up and not down.	Hand control defective.	Do not use the ceiling hoist and contact your Invacare provider for service.
	Strap length reference was lost.	Perform a system reset. Refer to 5.10 <i>Resetting the System, page 15.</i>
Lifting straps are twisted and do not straighten out.	Lifting hooks do not turn freely enough.	Check that the hook rotate freely and clean if necessary.
	Straps are frayed.	Do not use the ceiling hoist and contact your Invacare provider for service.
Lifting module emits excessive noise when activated.	Bearings, gearwheels or motor defective.	Do not use the ceiling hoist and contact your Invacare provider for service.
No indicator light on lifting module.	System does not have power or battery is discharged.	Charge the battery. Refer to 5.6 <i>Charging the Battery, page 14.</i>
The emergency strap does not turn the lift off	Ceiling hoist requires service.	Do not use the ceiling hoist and contact your Invacare provider for service.
The emergency lowering strap does not lower the lifting straps.	Ceiling hoist requires service.	Use the mechanical emergency lowering. Do not use the ceiling hoist and contact your Invacare provider for service.
Ceiling hoist does not traverse the rails.	Rails require service or cleaning.	Contact your Invacare provider for service.
	Trolley worn or damaged.	Trolley must be replaced. Contact your Invacare provider for service.
Ceiling hoist does not charge.	System has no power.	Ensure the hand control is correctly plugged to the charging station.  Leds blink to indicate charging is ongoing.
		Contact your Invacare provider for service.

10 Technical Data

10.1 Dimensions and weight

Dimensions

Robin® EVO	
Hoist length (A)	505 mm
Hoist width (B)	190 mm
Hoist height (C)	211 mm
Maximum lifting range (D)	2.5 m

Weights

Robin® EVO	
Maximum lifting capacity (Safe Working Load)	200 kg
Total weight without sling	9.5 kg

10.2 Electrical System

Robin® EVO	
Maximum current input	2.3 A
Protection class ¹	Entire device: IPx4 (ceiling hoist: IP24; hand control: IP44)
Insulation class	Class II equipment
	Type BF applied part Applied part complying with the specified requirements for protection against electrical shock according to IEC 60601-1.
Sound level	55 dB (A) [measured 1.0 m away from the device]
Number of lifts per charge	90 lifts of 0.5 m with 80 kg 60 lifts of 0.5 m with 200 kg
Intermittent (periodic motor operation)	10%, max. 2 min. / 18 min.
Battery	Li-Ion 36 V / 3.35 Ah
Charging time	Up to 7 hours at 20°C
Lifting speed — without load	4.2 cm/s
Lifting speed — with a 200 kg load	3.7 cm/s

Power Supply

Voltage output	36 V DC
Voltage supply	100 – 240 V AC, 50 /60 Hz, 1A
Protection class ¹	IP00

¹ See product label and label on each electric device for correct protection class. The lowest IP classification decides the overall classification of the device.

- IPx4: Protected against water splashed from any direction.
- IP24: Protected against objects larger than 12.5 mm and protected against water splashed from any direction.
- IP44: Protected against objects larger than 1 mm and protected against water splashed from any direction.

10.3 Environmental conditions

	Storage and transportation	Operation
Temperature	-10 °C to +50 °C	+5 °C to +40 °C
Relative humidity	20 % to 90 %	20 % to 90 % at 30 °C, not condensing
Atmospheric pressure	860 hPa to 1060 hPa	

Allow the product to reach operation temperature before usage:

- Warming up from minimum storage temperature might take up to 24 hours.
- Cooling down from maximum storage temperature might take up to 24 hours.

11 Electromagnetic compatibility (EMC)

11.1 General EMC information

Medical Electrical Equipment needs to be installed and used according to the EMC information in this manual.

This product has been tested and found to comply with EMC limits specified by IEC/EN 60601-1-2 for Class B equipment.

Portable and mobile RF communications equipment can affect the operation of this product.

Other devices may experience interference from even the low levels of electromagnetic emissions permitted by the above standard. To determine if the emission from this product is causing the interference, run and stop running this product. If the interference with the other device operation stops, then this product is causing the interference. In such rare cases, interference may be reduced or corrected by the following:

- Reposition, relocate, or increase the separation between the devices.

11.2 Electromagnetic emission

Guidance and manufacturer's declaration

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group I	This product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This product is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

11.3 Electromagnetic Immunity

Guidance and manufacturer's declaration

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Immunity test	Test / Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines; 100 kHz repetition frequency ± 1 kV for input / output lines; 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines	< 0% U_T for 0,5 cycle at 45° steps 0% U_T for 1 cycles 70% U_T for 25 / 30 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of this product requires continued operation during power mains interruptions, it is recommended that the product is powered from an uninterruptible power supply or a battery. U_T is the a. c. mains voltage prior to application of the test level.

Immunity test	Test / Compliance level	Electromagnetic environment – guidance
IEC 61000-4-11	< 5% U_T for 250 / 300 cycles	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 V 150 kHz to 80 MHz 6 V in ISM & amateur radio bands	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this product is used exceeds the applicable RF compliance level above, this product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating this product.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz 385 MHz - 5785 MHz test specifications for immunity to RF wireless communication equipment refer to table 9 of IEC 60601-1-2	Interference may occur in the vicinity of equipment marked with the following symbol:  Portable and mobile RF communications equipment should be used no closer than 30 cm to any part of this product including cables.

 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

11.4 EMC test specifications

IEC 60601-1-2 — Table 9

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	28
710 745 788	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	28
1720 1845 1970	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	28
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	28
5240 5500 5785	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	9

 If necessary to achieve the immunity test level, the distance between the transmitting antenna and the Medical Electrical Equipment or System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal.

^{c)} As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz.

While it does not represent actual modulation, it would be worst case.

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Invacare AB
Invacare UK Operations Limited